



## Biomol BLUEplus prestained Protein Ladder (10-180 kDa)

**Catalog No:** 10008  
**Lot No:** XXXXX  
**Supplied as:** 500 µl  
**Stability:** store at 4°C for 3 months or at -20°C for 24 months

### Background

The Biomol BLUEplus prestained Protein Ladder is a three-color protein standard with 10 pre-stained proteins covering a wide range molecular weights from 10 to 180 kDa. Proteins are covalently coupled with a blue chromophore except for two reference bands (one green and one red band at 25 kDa and 75 kDa respectively) when separated on SDS-PAGE (Tris-glycine buffer). The Biomol BLUEplus prestained Protein Ladder is designed for monitoring protein separated during SDS-polyacrylamide gel electrophoresis, verification of Western transfer efficiency on membranes (PVDF, nylon, or nitrocellulose) and for approximate sizing of proteins. The ladder is supplied in gel loading buffer and is ready to use. Do not heat, dilute, add reducing agent before loading.

### Formulation


Approximately 0.2 – 0.4 mg/ml of each protein in buffer (20 mM Trisphosphate pH 7.5 at 25°C), 2% SDS, 0.2 mM Dithiothreitol, 3.6 M Urea, and 15% (v/v) Glycerol.

### Protocol for use in SDS-PAGE

- Thaw the ladder either at room temperature or at 37-40°C for a few minutes to dissolve precipitated solids. Do not boil.
- Mix thoroughly to ensure the solution is homogeneous.
- Load the following volumes of the ladder on SDS-PAGE:

**Mini gel:** 3 µl per well for SDS-PAGE, 1.5 µl for WB

**Large gel:** 6 µl per well for SDS-PAGE, 3 µl for WB



### Usage

This product is offered by Biomol for research purposes only. Not for diagnostic purposes or human use. It may not be resold or used to manufacture commercial products without written approval of Biomol GmbH.

### CONTACT US TODAY

BIOMOL GmbH • Kieler Straße 303a • 22525 Hamburg • Germany • [info@biomol.de](mailto:info@biomol.de) • [www.biomol.de](http://www.biomol.de)

Fon: +49 (0)40-853 260 0 • TOLL FREE IN GERMANY: Fon: 0800-246 66 51