



Leptin Antagonist Triple Mutant, mouse recombinant (rmLeptin-tm)

Catalog No: 86418
Lot No: XXXXX
Source: *E. coli*
Synonyms:

Description

Leptin antagonist triple mutant mouse recombinant is a single non-glycosylated polypeptide chain containing 146 amino acids and additional Ala at N-terminus and having a molecular mass of approx. 16 kDa. Mouse leptin was mutated, resulting in L39A/D40A/F41A mutant. Leptin antagonist triple mutant was purified by proprietary chromatographic techniques.

Physical Appearance

White lyophilized (freeze-dried) powder.

Formulation

The protein was lyophilized from a concentrated (0.65 mg/ml) solution with 0.003 mM NaHCO₃.

Solubility

It is recommended to reconstitute the lyophilized leptin antagonist triple mutant in sterile water or sterile 0.4% NaHCO₃ adjusted to pH 8-9, not less than 100 µg/ml, which can then be further diluted with other aqueous solutions.

Stability

Lyophilized leptin antagonist triple mutant, although stable at room temperature for several weeks, should be stored desiccated below -18°C. Upon reconstitution at > 0.1 mg/ml leptin mutant and up to 2 mM and filter sterilization LEP mutant can be stored at 4°C or even room temperature for several weeks making it suitable for long term infusion studies using osmotic pumps. At lower concentration addition of a carrier protein (0.1% HSA or BSA) is suggested. Please prevent freeze-thaw cycles.

Purity

Greater than 99.0% as determined by (a) Gel filtration analysis, (b) Analysis by SDS-PAGE.

Amino Acid Sequence

The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Val-Pro-Ile-Gln.

Activity

Leptin antagonist triple mutant mouse recombinant is capable of inhibiting leptin-induced proliferation of BAF/3 cells stably transfected with the long form of human leptin receptor. It also inhibits various leptin effects in several in vitro bioassays.

Usage

This product is offered by Biomol for research purposes only. Not for diagnostic purposes or human use. It may not be resold or used to manufacture commercial products without written approval of Biomol GmbH.

CONTACT US TODAY

BIOMOL GmbH • Kieler Straße 303a • 22525 Hamburg • Germany • info@biomol.de • www.biomol.de

Fon: +49 (0)40-853 260 0 • TOLL FREE IN GERMANY: Fon: 0800-246 66 51