

YOU HAVE THE VISION, WE HAVE THE SUBSTANCE.



# Flt3-Ligand, human recombinant (rHuFlt3-Sf9)

Catalog No:87412Lot No:XXXXSource:Baculovirus Sf9 cellsSynonyms:Fms-related tyrosine kinase 3 ligand, FLK2, STK1, CD135, Stem Cell Tyrosine Kinase 1, FLT3LG, Flt3

# Background

FLT3 ligand is a receptor for the fl cytokine has a tyrosine-protein kinase activity & a growth factor that regulates proliferation of early hematopoietic cells. Flt3-Ligand synergizes with other CSFs and interleukins to induce growth and differentiation.

## Description

Flt3-Ligand human recombinant produced in Sf9 insect cells is a glycosylated, polypeptide chain containing 155 amino acids and is migrating on SDS-PAGE at 19 kDa. Flt3-Ligand is purified by proprietary chromatographic techniques.

#### **Physical Appearance**

Sterile filtered white lyophilized (freeze-dried) powder.

#### Formulation

The protein was lyophilized with no additives.

## Solubility

It is recommended to reconstitute the lyophilized Flt3-Ligand in sterile 18 M $\Omega$ -cm H<sub>2</sub>O not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.

#### Stability

Lyophilized Flt3-Ligand, although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Flt3-L should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

#### Purity

Greater than 95.0% as determined by (a) Analysis by RP-HPLC, (b) Analysis by SDS-PAGE.

#### **Amino Acid Sequence**

The sequence of the first five N-terminal amino acids was determined and was found to be Thr-Gln-Asp-Cys-Ser.

#### Activity

The ED50, calculated by the dose-dependent stimulation of the proliferation of human MUTZ-2 cells is <3.0 ng/ml.

## Usage

This product is offered by Biomol for research purposes only. Not for diagnostic purposes or human use. It may not be resold or used to manufacture commercial products without written approval of Biomol GmbH.

# **CONTACT US** TODAY

BIOMOL GmbH • Kieler Straße 303a • 22525 Hamburg • Germany • info@biomol.de • www.biomol.de Fon: +49 (0)40-853 260 0 • TOLL FREE IN GERMANY: Fon: 0800-246 66 51