Cat IgG-heavy and light chain Antibody

Cat IgG-heavy and light chain Antibody			
Rabbit Polyclonal Antiserum			
Catalog No. A20–115			BETHYL
Lot No. A20–115–2			LABORATORIES, INC
APPLICATIONS	IEP, DD		
SPECIES REACTIVITY	Cat		
ISOTYPE	IgG		
AMOUNT	2 ml		
STORAGE/SHELF LIFE	2 - 8° C / 2 years from date of receipt		
PHYSICAL STATE Liquid			
BUFFER	Serum containing 0.09% Sodium Azide		
ORIGIN	USA		
PRODUCTION PROCEDURES	By immunoelectrophoresis, the antiserum reacts specifically with Fab fragment of cat IgG and with light chains common to other cat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.		
APPLICATIONS	Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.		
	IEP	Neat	
	DD	Neat	
APPLICATION NOTES	For use in precipitin gel Ouchterlony.	reactions; such as, immunoelectrophoresis and	double diffusion
ADDITIONAL INFO	https://www.bethyl.com/product/A20-115 Use the link above to view SDS, a current list of citations, and other product specific information.		

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc. Eric McIntush, PhD | Chief Scientific Officer Date: December 3, 2018

For in vitro laboratory use only. Not for any clinical, therapeutic or diagnostic use in humans or animals. Not for human or animal consumption. This product may not be resold or modified for resale without the prior written approval of Bethyl Laboratories, Inc. The information provided in this data sheet is believed to be correct but does not purport to be all-inclusive and is intended to be used as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state or local laws or regulations.