





#### Detection of human CtIP by immunoprecipitation.

*Samples:* Whole cell lysate (1 mg for IP; 1/4 of IP loaded) from MCF-7 cells. *Antibodies:* Affinity purified rabbit anti-CtIP antibody BL1913 (Cat. No. A300-487A) used at 3  $\mu\text{g}/\text{mg}$  lysate for IP. CtIP was also immunoprecipitated with rabbit anti-CtIP antibody BL1914 (Cat. No. A300-488A) using 3  $\mu\text{g}/\text{mg}$  lysate. immunoprecipitatesd CtIP was blotted using BL1914 at 1  $\mu\text{g}/\text{ml}$ . *Detection:* Chemiluminescence with an exposure time of 30 seconds.