



Detection of human BEGAIN by western blot of immunoprecipitates. *Samples:* Whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells. *Antibodies:* Affinity purified rabbit anti-BEGAIN antibody A302-792A used for IP at 10 $\mu\text{g}/\text{mg}$ lysate. BEGAIN was also immunoprecipitated by rabbit anti-BEGAIN antibody A302-793A, which recognizes a downstream epitope. For blotting immunoprecipitated BEGAIN, A302-793A was used at 1 $\mu\text{g}/\text{ml}$. *Detection:* Chemiluminescence with an exposure time of 30 seconds.