





**Detection of human Kaiso by western blot and immunoprecipitation.**

*Samples:* Whole cell lysate from HEK293T (15 and 50  $\mu$ g for WB; 1 mg for IP, 20% of IP loaded), Jurkat (J; 50  $\mu$ g) and HeLa (H; 50  $\mu$ g) cells.

*Antibodies:* Affinity purified rabbit anti-Kaiso antibody A303-557A used for WB at 0.04  $\mu$ g/ml (A) and 1  $\mu$ g/ml (B) and used for IP at 6  $\mu$ g/mg lysate. Kaiso was also immunoprecipitated by rabbit anti-Kaiso antibody A303-558A, which recognizes a downstream epitope.

*Detection:* Chemiluminescence with exposure times of 3 minutes (A) and 30 seconds (B).