





**Detection of human PFKFB3 by western blot of immunoprecipitates.** *Samples:* Whole cell lysate (0.5 to 1.0 mg per IP reaction; 20% of IP loaded) from Jurkat cells. *Antibodies:* Affinity purified rabbit anti-PFKFB3 antibody A304-248A (lot A304-248A-1) used for IP at 6 µg per reaction. PFKFB3 was also immunoprecipitated by rabbit anti-PFKFB3 antibody A304-249A. For blotting immunoprecipitated PFKFB3, A304-249A was used at 1 µg/ml. *Detection:* Chemiluminescence with an exposure time of 30 seconds.