





**Detection of human TPP2 by western blot of immunoprecipitates.** *Samples:* Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from 293T cells prepared using NETN lysis buffer. *Antibodies:* Affinity purified rabbit anti-TPP2 antibody A304-758A (lot A304-758A-1) used for IP at 6 µg per reaction. TPP2 was also immunoprecipitated by rabbit anti-TPP2 antibody A304-757A. For blotting immunoprecipitated TPP2, A304-757A was used at 0.4 µg/ml. *Detection:* Chemiluminescence with an exposure time of 10 seconds.