





**Detection of human EXOC4/SEC8 by western blot of immunoprecipitates.** *Samples:* Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from Jurkat cells prepared using NETN lysis buffer. *Antibody:* Affinity purified rabbit anti-EXOC4/SEC8 antibody A304-785A (lot A304-785A-1) used for IP at 6 µg per reaction. EXOC4/SEC8 was also immunoprecipitated by rabbit anti-EXOC4/SEC8 antibody A304-786A. For blotting immunoprecipitated EXOC4/SEC8, A304-786A was used at 1 µg/ml. *Detection:* Chemiluminescence with an exposure time of 30 seconds.