



Detection of human NUDCD3 by western blot of immunoprecipitates. *Samples:* Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from 293T cells prepared using NETN lysis buffer. *Antibody:* Affinity purified rabbit anti-NUDCD3 antibody A304-865A (lot A304-865A-1) used for IP at 6 µg per reaction. NUDCD3 was also immunoprecipitated by rabbit anti-NUDCD3 antibody A304-866A. For blotting immunoprecipitated NUDCD3, A304-865A was used at 0.4 µg/ml. *Detection:* Chemiluminescence with an exposure time of 10 seconds.