





**Detection of human PEX1 by western blot of immunoprecipitates.** *Samples:* Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from Jurkat cells prepared using NETN lysis buffer. *Antibodies:* Affinity purified rabbit anti-PEX1 antibody A305-386A (lot A305-386A-1) used for IP at 6 µg per reaction. PEX1 was also immunoprecipitated by rabbit anti-PEX1 antibody A305-387A. For blotting immunoprecipitated PEX1, A305-386A was used at 1 µg/ml. *Detection:* Chemiluminescence with an exposure time of 3 minutes.