



Detection of human MTMR6 by western blot of immunoprecipitates. *Samples:* Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. *Antibodies:* Affinity purified rabbit anti-MTMR6 antibody A305-496A (lot A305-496A-1) used for IP at 6 µg per reaction. MTMR6 was also immunoprecipitated by rabbit anti-MTMR6 antibody BL21798. For blotting immunoprecipitated MTMR6, A305-496A was used at 0.4 µg/ml. *Detection:* Chemiluminescence with an exposure time of 10 seconds.