



Custom Mouse Primer Library

Catalog No: CMPL-1
Supplier: RealTimePrimers
Lot No: XXXXX
Supplied as: solid
Stability: store at -20°C

Description

Primers are provided in a dry format and may be reconstituted with 40ul of dH₂O to obtain a 50 uM solution. Enough for 1000 PCR arrays based on a 20ul reaction at a final conc of 0.1 uM. Choose up to 96 genes per plate.

Usage

This product is offered by Biomol for research purposes only. Not for diagnostic purposes or human use. It may not be resold or used to manufacture commercial products without written approval of Biomol GmbH.