



Diacclone

DIACLONE Immunology Products

IL-7 ELISA Kit

Human IL-7 ELISA Kit

CATALOGUE N° :	850.040.048, 850.040.096, 850.040.192
SPECIFICITY :	Recognizes both natural and recombinant human IL-7
RANGE :	6.25 pg / ml - 200 pg / ml
SENSITIVITY :	< 3 pg / ml
INCUBATION :	2h40
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human cytokines
KIT SIZE :	48 (Discovery), 96, 192 tests
REFERENCES :	Bikker, A. et al., Arthritis Rheum.,2010; 62(4): 969-77 Butler, M. O. et al., PLoS One,2012; 7(1): e30229 Churchman, S. M. and F. Ponchel, Rheumatology, 2008;47:753-759 Kang, I. et al., J. Immunol., 2004; 173 : 673 - 681 Levy, Y. et al., J Clin Invest.,2009 ; 119(4): 997-1007 Perales, M-A et al., Blood,2012;120:4882-4891 Plana, M. et al., Immunology,2011; 133(3): 318-28 Sereti, I. et al., Blood,2009;113(25):6304-6314 Sportes, C. et al., Clin. Cancer Res.,2010;16(2): 727-735 Thiebaut, R. et al., Clinical Infectious Diseases,2016;62(9): 1178-1185 Thiebaut, R. et al., Clinical Infectious Diseases,2016;62(9): 1178-1185 van Roon, J. A. et al., Ann Rheum Dis.,2007; 66(5): 664-9. van Roon, J. A. et al., Ann Rheum Dis.,2003; 62(2): 113-9.
KIT CONTENT :	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Calibrated to WHO 90/530

Version 11 - 12.16

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985

25020 Besançon cedex - France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

info@diacclone.com



www.diacclone.com