



Diacclone

DIACLONE Immunology Products

IL-13 ELISA Kit

Human IL-13 ELISA kit

CATALOGUE N° :	850.080.048, 850.080.096, 850.080.192
SPECIFICITY :	Recognizes both natural and recombinant human IL-13
RANGE :	3.12 pg / ml - 100 pg / ml
SENSITIVITY :	< 1.5 pg / ml
INCUBATION :	2h10
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human cytokines
KIT SIZE :	48 (Discovery), 96, 192 tests
REFERENCES :	Abbal C. et al., Int. Immunol., 1999; 11(9) : 1451 - 1462 Antczak, A. et al., International Journal of Immunopathology and Pharmacology, 2016: 0394632015623794 Cameron, L. et al., J Immunol., 2006; 177(12): 8633-42. de Nadai, P. et al., J Immunol., 2006; 176(10): 6286-93. Fahy O. et al., J. Immunol., 2002; 168(11) : 5912 - 5919 Finan, C. et al., PLoS One, 2008; 3(10): e3485 Gafvelin, G. et al., J Biol Chem., 2007; 282(6): 3778-87. Halonen, M. et al., J. Immunol., 2009; 182(5): 3285-3293. Ismaili, J. et al., Clin Exp Immunol., 2003 ; 133(3): 414-21 Leung S. et al., J. Biomol Screen, 2005; 10(2): 157 - 167 Pene, J. et al., J. Immunol., 2008; 180(11): 7423-7430. Perrin, Y. et al., Clin Transl Allergy, 2014; 4(1): 1 Piccoli, L. et al., Clin Dev Immunol., 2012: 483935 Pitta, M. G. et al., J Clin Invest., 2009; 119(8): 2379-87 Saverino D. et al., J. Immunol., 2002; 168(1): 207 - 215 Schmitt, C. et al., J Leukoc Biol., 2000; 68(6): 836-44. Tsvetkova-Vicheva, V. M. et al., Clin Transl Allergy, 2014; 4(1): 3 Vekemans, J. et al., Clin Exp Immunol., 2002 ; 127(3): 495-8 Woerty, G. et al., J Leukoc Biol., 2002; 72(4): 769-79.
KIT CONTENT :	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Calibrated to WHO 94/622

Version 10 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com