



Diacclone


DIACLONE Immunology Products

CD130 ELISA Kit

Human CD130 ELISA kit

CATALOGUE N° :	850.520.048, 850.520.096, 850.520.192
SPECIFICITY :	Recognizes both natural and recombinant human soluble gp130
RANGE :	56.25 pg / ml - 1800 pg / ml
SENSITIVITY :	< 15 pg / ml
INCUBATION :	1h40
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human soluble receptors No interference with neither soluble gp80 nor IL-6
KIT SIZE :	48 (Discovery), 96, 192 tests

REFERENCES :	Arranz-Valsero, I. et al., Mol Vis., 2013; 19: 85-99 Corbi P. et al., Eur J Cardiothorac Surg., 2000; 18: 98 - 103 Dore, P. et al., Clin Exp Immunol., 1997; 107(1): 182-8 Gougelet, A. et al., PLoS One, 2009; 4(11): e8026 Muller-Newen, G. et al., J Immunol., 1998; 161 (11): 6347-55 Santer, F. R. et al., Endocr. Relat. Cancer, 2010; 17(1): 241-253
---------------------	--



KIT CONTENT :	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.
----------------------	--

Version 7 - 12.16

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com