



Diacclone

DIACLONE Immunology Products

CD62E ELISA Kit

Human CD62E ELISA kit

CATALOGUE N° :	850.530.048, 850.530.096, 850.530.192
SPECIFICITY :	Recognizes both natural and recombinant human soluble ELAM-1
RANGE :	1.0 ng / ml - 32 ng / ml
SENSITIVITY :	< 0.5 ng / ml
INCUBATION :	2h10
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human adhesion molecules
KIT SIZE :	48 (Discovery), 96, 192 tests

REFERENCES :	Adepu, S. et al., Am J Physiol Renal Physiol.,2015; 309(2): F137-145 Altun, I. et al., J Clin Med Res.,2014; 6(5): 354-61 Bian, C. et al., J Zhejiang Univ Sci B.,2005; 6(7): 631-6 Branger, J. et al., Blood,2003; 101(11): 4446-8. Chirico, E.N. et al., J Appl Physiol, 2012; 112: 1445 - 1453 Domienik-Karłowicz, J. et al.,Obes Surg.,2014; 25(3): 486-93 Engelen, L., et al., Eur. J. Endocrinol.,2011; 164(3): 371-379 Gombocz, K. et al., Crit Care,2007; 11(4): R87 Hartweg, J. et al., Clin. Chem.,2007; 53(10): 1858-1860. Kuenen B.C. et al., Clin. Cancer Res., 2005; 11(17): 6240 - 6246 Larsson, S. et al., Int J Chron Obstruct Pulmon Dis.,2011 ; 6: 83-8 le Poole, C. Y. et al.,Perit. Dial. Int.,2013;33(2): 208-212 Lund, S.S. Et al., Eur. J. Endocrinol.,2008;158(5): 631-641. Manneras-Holm, L. et al., J. Clin. Endocrinol. Metab.,2011; 96(4): 1068-1076. Parchure, N. et al., Circulation,2002; 105(11): 1298-303. Pelliccione, F. et al., Asian J Androl.,2014;16(2):290-4 Persson, F. et al., Diabetes,2006; 55(12): 3550-5. Renckens R. et al., Arterioscler. Thromb. Vasc. Biol., 2004; 24(3): 483 - 488 Trappenburg, M. C. et al., Haematologica,2009; 94(7): 911-918. Van Dam, B. et al.,Cardiovasc Res.,2003; 57(2): 563-571. Van Ree R M. et al., Nephrol. Dial. Transplant., 2005; 20(11):2524-31 Verier, C. M. P. et al., J. Nutr.,2011; 141(3): 411-417 Vincent, L. et al., Am J Physiol Heart Circ Physiol.,2009: 00812.2009. Yao, G. H. et al., Lupus,2008;17(8): 720-726.
---------------------	--

KIT CONTENT :	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.
----------------------	--

Version 13 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com