



Diacclone

## DIACLONE Immunology Products

### CD62L ELISA Kit

Human CD62L ELISA kit

<b>CATALOGUE N° :</b>	850.650.048, 850.650.096, 850.650.192
<b>SPECIFICITY :</b>	Recognizes both natural and recombinant human soluble LAM-1
<b>RANGE :</b>	1.6 ng / ml - 60 ng / ml
<b>SENSITIVITY :</b>	< 1 ng / ml
<b>INCUBATION :</b>	1h45
<b>SAMPLE TYPES :</b>	Serum Plasma Cell culture supernatant
<b>SAMPLE SIZE :</b>	100 µl
<b>CROSS REACTION :</b>	No cross reactivity with other human adhesion molecules
<b>KIT SIZE :</b>	48 (Discovery), 96, 192 tests  Capo, X. et al., <i>Nutrients</i> , 2016; 8(10). Fabre, S. et al., <i>J. Immunol.</i> , 2008; 181(5): 2980-2989. Nebor, D. et al., <i>Haematologica</i> , 2011; 96(11): 1589-1594 Reiss, K. et al., <i>J. Biol. Chem.</i> , 2011; 286: 26931 - 26942 Shpacovitch V.M et al., <i>J. Leukoc. Biol.</i> , 2004; 76(2): 388 - 398 Tripette, J. et al., <i>Br. J. Sports Med.</i> , 2010; 44(4): 232-237 Verier, C. M. P. et al., <i>J. Nutr.</i> , 2011; 141(3): 411-417
<b>REFERENCES :</b>	
<b>KIT CONTENT :</b>	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Version 8 - 06.17

For research use only

### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

#### DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod  
- BP 1985  
25020 Besançon cedex - France  
Tel: +33 (0)3 81 41 38 38  
Fax: +33 (0)3 81 41 36 36  
info@diacclone.com



[www.diacclone.com](http://www.diacclone.com)