



Diacclone

DIACLONE Immunology Products

TNF- α ELISA Set

Human TNF- α ELISA Set

CATALOGUE N° :	851.570.001, 851.570.005, 851.570.010, 851.570.015, 851.570.020
SPECIFICITY :	Recognizes both natural and recombinant human TNF- α
RANGE :	Recommended from 25 pg / ml to 800 pg / ml
SENSITIVITY :	< 10 pg / ml
INCUBATION :	From sample to end 2h45
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 μ l
CROSS REACTION :	No cross reactivity with other human cytokines Cross reactivity with macac TNF- α No interference with soluble TNFR I and TNFR II
KIT SIZE :	1X96 (Discovery), 5X96, 10X96, 15X96, 20X96 tests Filaci, G. et al., J. Immunol.,2007; 179(7): 4323-4334 Fomander, L. et al.,PLoS One,2013; 8(12): e83089 Haller D. et al., Clin. Diagn. Lab. Immunol., 2002; 9(3): 649 - 657 Harel, E. et al., PLoS One,2011; 6(9): e24202
REFERENCES :	Noizat-Pirenne, F. et al., Haematologica,2007; 92(12): e132-135. Popov, A. et al., J Clin Invest,2006; 116(12): 3160-70. Popov, A. et al.,J. Immunol.,2008; 181(7): 4976-4988. Sabour Alaoui, S. et al.,PLoS One,2012;7(3): e33609 Zhang, Q. et al., Infect Immun.,2006; 74(8): 4735-43.
KIT CONTENT :	Capture antibody,Biotinylated detection antibody,Standard,Streptavidin-HRP,TMB

Version 8 - 06.17 For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com