



Diacclone

DIACLONE Immunology Products

Anti-Human CD41a Unconjugated

Clone: HIP8

CATALOGUE N° :	853.223.020
CLONE :	HIP8
ISOTYPE :	Mouse IgG1
SPECIFICITY :	The antibody HIP8 reacts with CD41a antigen. CD41 (GPIIb) is a transmembrane glycoprotein (integrin family) composed of two chains GPIIb α and GPIIb β .
WORKSHOP :	IV
QUANTITY :	200 μ g / 200 μ l
FORMAT :	The monoclonal is presented in phosphate buffered saline (PBS) with 15 mM sodium azide, approx. pH 7.4
WORKING DILUTION :	We recommend that each laboratory determine an optimum working titre for use in its particular application.
STORAGE :	Store at +2-8°C. For longer storage, freeze aliquots
APPLICATION :	Flow Cytometry. Immunohistochemistry (frozen sections).
OTHER :	CD41 is non-covalently associated with CD61 (GPIIIa) to form the GPIIb/IIIa (CD41/CD61) complex
REFERENCES :	Bao CX, et al, Chinese J. of Hematology., 8(2):65 Bao CX, et al, Chinese J. of Hematology., 13(2):66 Chen Z., et al, Chinese Science Bulletin., 24:1902 Wu XW., et al., National Medical J. of China., 1989; 69(8):427 Leucocyte Typing IV. Knapp W et al. (Eds.), Oxford University Press (1989). Leucocyte Typing V. Schlossman S. et al. (Eds.), Oxford University Press (1995). Leukocyte Typing VI. Kishimoto T. et al. (Eds.), Garland Publishing Inc. (1997).

☛ Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

Version 3 – 03.16

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
– BP 1985

25020 Besançon cedex – France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

info@diacclone.com



www.diacclone.com