


## Anti-Human Galectin-8 Unconjugated

### PRODUCT SPECIFICATIONS

Catalogue N°	857.013.020 - 200 tests / 2ml
Target species	Human
Specificity	Recognizes both natural and recombinant human galectin-8
Clone	B-Z18
Application	Flow Cytometry
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells
Immunsation	Recombinant human Galectin-8
Quantity	200 Tests / 2ml (Discovery Size also available please enquire)
Isotype	Mouse IgG2a Kappa light chain
Format	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Working Dilution	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
Synonym	LGALS8 Gal-8 PCTA-1



A typical staining pattern of YT (T cells) with the B-Z18 monoclonal antibody  
Black peak: Isotype control  
Red peak: Labelled cells

#### Discover Diaclone

With over 30 years experience and extensive expertise, we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diaclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity and quality of our product range makes Diaclone a clear choice to :

#### Fast Track Your Research

##### Diaclone products



6 rue Dr Jean-François-Xavier Girod  
BP 1985  
25020 Besançon cedex  
France

Tel : +33 (0)3 81 41 38 38  
Fax : +33 (0)3 81 41 36 36

Email : [Info@diaclone.com](mailto:Info@diaclone.com)  
[www.diaclone.com](http://www.diaclone.com)

Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

