



Diacclone

DIACLONE Immunology Products

Anti-Human CD81 Azide Free

Clone: TS81

CATALOGUE N° :	857.780.000
CLONE :	TS81
ISOTYPE :	IgG2a
IMMUNISATION :	Jurkat cell line
HYBRIDOMA :	Myeloma X63/AG8. x Balb/c spleen cells
SPECIFICITY :	Recognises the TAPA-1 antigen, a 23 kDa (smear) protein
QUANTITY :	200 µg / 200 µl
FORMAT :	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
STORAGE :	Stable at +2-8°C for 12 months For longer storage, freeze aliquots
APPLICATION :	Suitable for flow cytometry, immunoprecipitation and Western-Blot (under non-reducing conditions)

REFERENCES :	Abache, T. et al., J Cell Biochem, 2007; 102(3): 650-64 Charrin S. et al., Eur. J. Immunol., 2003; 33: 2479-2489 Charrin S. et al., Biochem. J., 2003; 373: 409-421 Charrin S. et al., Biochem. Biophys. Res. Commun., 2003; 304: 107-112 Charrin S. et al., FEBS Lett., 2002; 516: 139-144 Charrin S. et al., J. Biol. Chem., 2001; 276: 14329-14337 Serman, A. et al., Nucleic Acids Res., 2007; 35 (14): 4715-4727. Silvie, O. et al., J Cell Sci., 2006; 119(Pt 10): 1992-2002 Verweij, F. J. et al., Embo J., 2011; 30(11): 2115-29 Yalaoui, S. et al., PLoS Pathog., 2008; 4(2): e1000010
---------------------	--

97.4

58.1

39.8

29

20



1,3 105 Raji cells were lysed (lysis buffer containing 1% Brij97 or triton)

Migration SDS/page 5-15% gel under non-reducing conditions

Version 7 - 12.16

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985

25020 Besançon cedex - France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

info@diacclone.com



www.diacclone.com