



Diacclone

DIACLONE Immunology Products


Anti-Human CD82 Azide Free

Clone: TS82b

CATALOGUE N° :	857.790.000
CLONE :	TS82b
ISOTYPE :	Mouse IgG1
IMMUNISATION :	Jurkat and HEL cell lines
HYBRIDOMA :	Myeloma X63/AG8. x Balb/c spleen cells
SPECIFICITY :	Recognises the CD82 antigen, a 60 kDa protein
QUANTITY :	200 µg / 200 µl
FORMAT :	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
STORAGE :	Stable at +2-8°C for 12 months DO NOT FREEZE.
APPLICATION :	Suitable for flow cytometry, immunoprecipitation and Western-Blot (under non-reducing conditions)

REFERENCES :

Bari, R. et al., Am. J. Pathol., 2009; 174(2): 647-660.
 Charrin S. et al., Biochem. J., 2003; 373(Pt 2): 409-421
 Charrin S. et al., Eur. J. Immunol., 2003; 33(9): 2479-2489
 Charrin S. et al., FEBS Lett., 2002; 516(1-3): 139-144
 Charrin S. et al., J. Biol. Chem., 2001; 276(17): 14329-14337
 Charrin, S. et al., Biochem Biophys Res Commun., 2003; 304(1): 107-12.
 Gellersen, B. et al., Am. J. Pathol., 2007 ;170(1):126-139.
 Gellersen, B. et al., Hum. Reprod., 2010; 25(4): 862-873
 Gellersen, B., et al., Biol Reprod., 2013; 88(3):80
 Gonzalez, M. et al., Mol. Hum. Reprod., 2011; 17(7): 421-433
 Poloso, N. J. et al., J Immunol., 2006; 177(8): 5451-8.
 Ruseva, Z. et al., Exp Cell Res., 2009; 315(10): 1759-71.
 Samalecos, A et al., Reprod Biol Endocrinol., 2009; 7: 76.
 Xu, C. et al., Faseb J., 2009; fj.08-123414



1,3 105 HEL cells were lysed (lysis buffer containing 1% Brij97 or triton)
 Migration SDS/page 5-15% gel under non-reducing conditions

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS

-  Monoclonal Antibodies
-  ELISA
-  ELISpot
-  DIAPlex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
 – BP 1985
 25020 Besançon cedex – France
 Tel: +33 (0)3 81 41 38 38
 Fax: +33 (0)3 81 41 36 36
 info@diacclone.com



www.diacclone.com