

Murine IL-2 ELISpot Kit

PRODUCT SPECIFICATIONS

Catalogue N°	862.001.001PC - 1 x 96 Discovery (pre-coated plate) 862.001.002PC - 2 x 96 (pre-coated plates) 862.001.005PC - 5 x 96 (pre-coated plates)
Target species	Murine
Specificity	Recognizes natural murine IL-2
Incubation	3h after cell stimulation
Cross Reaction	No cross reactivity with other murine cytokines
Kit Content	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.

REFERENCES

Pollard, C. et al., Mol Ther.,2013; 21(1): 251-9 - [Pubmed link](#) 

BACKGROUND

IL-2 is a powerful immunoregulatory lymphokine produced by T-Cells in response to antigenic or mitogenic activation. IL-2 stimulates growth and differentiation of B-Cells, most T-cells, NK cells, monocytes and macrophages.

Mature IL-2 is a 15.4kDa globular glycoprotein containing 133 amino acid residues including one intrachain disulfide bond between residues 58 and 105.

Apart from its most important role to mediate antigen-specific T-lymphocyte proliferation, IL-2 modulates the expression of interferon-g and major histocompatibility antigens.

Alterations in the ability of T-cell to synthesize IL-2 have been observed in physiologic and pathologic states. Currently, IL-2 is used to enhance the immune system of patients for the treatment of cancer and infectious disease.

Version 6 - 01.21

For research use only

Discover Diaclone

With over 30 years experience and extensive expertise, we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diaclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity and quality of our product range makes Diaclone a clear choice to :

Fast Track Your Research

Diaclone products

 **Monoclonal Antibodies**

 **ELISA**

 **ELISpot**

 **Multiplex**

 **Diaclone SAS**

6 rue Dr Jean-François-Xavier Girod
BP 1985
25020 Besançon cedex
France

Tel : +33 (0)3 81 41 38 38
Fax : +33 (0)3 81 41 36 36

Email : Info@diaclone.com

www.diaclone.com