



Diacclone

## DIACLONE Immunology Products

### IFN- $\gamma$ ELISpot Kit

Murine IFN- $\gamma$  ELISpot Enzymatic : Precoated PVDF plates

|                         |   |
|-------------------------|---|
| <b>CATALOGUE N° :</b>   | 862.031.001PC, 862.031.005PC  |
| <b>SPECIFICITY :</b>    | Recognizes natural murine IFN- $\gamma$   |
| <b>INCUBATION :</b>     | 3 h after cell stimulation  |
| <b>CROSS REACTION :</b> | No cross reactivity with other murine cytokines   |
| <b>KIT SIZE :</b>       | 1x96 (Discovery), 5x96 tests  |
| <b>REFERENCES :</b>     | Musson, J. A. et al., Infect. Immun., 2010; IAI.00195-10<br>Pollard, C. et al., Mbl Ther., 2013; 21(1): 251-9<br>Ramakrishna L. et al., J. Virol., 2004; 78(17): 9174 - 9189<br>Someya K. et al., J. Virol., 2004; 78(18): 9842 - 9853<br>Tsunoda, I. et al., J Virol., 2005; 79(23): 14640-6.<br>Yan, L. et al., Infect. Immun., 2008; IAI.00245-08. |
| <b>KIT CONTENT :</b>    | Precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.   |

Version 6 - 12.16

For research use only

### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

#### DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod  
- BP 1985

25020 Besançon cedex - France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

info@diacclone.com



[www.diacclone.com](http://www.diacclone.com)