



Diacclone

## DIACLONE Immunology Products

### Anti-Human CD126 Capture Antibody

Clone: B-N12

<b>CATALOGUE N° :</b>	879.510.001
<b>CLONE :</b>	B-N12
<b>ISOTYPE :</b>	Mouse IgG1
<b>IMMUNISATION :</b>	gp80 transfected CHO cell line
<b>HYBRIDOMA :</b>	Myeloma X63/AG.8653 x Balb/c spleen cells
<b>SPECIFICITY :</b>	Recognises the Interleukin-6 Receptor alpha subunit, IL-6Ra, IL-6R, a 80 kDa protein
<b>QUANTITY :</b>	1.0 mg / 1.0 ml
<b>FORMAT :</b>	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
<b>STORAGE :</b>	Stable at +2-8°C for 12 months For longer storage, freeze aliquots
<b>APPLICATION :</b>	<b>ELISA Capture Antibody</b> This Antibody can be used as Capture in a human CD126 sandwich Immunoassay to detect human CD126 in combination with biotinylated human CD126 Detection Antibody (Cat n° 879.510.002). The suggested coating concentration range below should be optimised by each laboratory for each application. <b>ELISA 0,3-1µg/ml</b> General ELISA protocol is available at <a href="http://www.diacclone.com">www.diacclone.com</a>

Version 2 - 03.16

For research use only

#### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

#### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

#### DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod  
– BP 1985

25020 Besançon cedex – France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

[info@diacclone.com](mailto:info@diacclone.com)



[www.diacclone.com](http://www.diacclone.com)