



Diacclone

## DIACLONE Immunology Products

### Anti-Human CD141 Detection Antibody

Clone: B-S24

<b>CATALOGUE N° :</b>	879.720.002
<b>CLONE :</b>	B-S24
<b>ISOTYPE :</b>	Mouse IgG1
<b>IMMUNISATION :</b>	Eahy 926 cell line
<b>HYBRIDOMA :</b>	Myeloma XG3/AG.8653 x Balb/c spleen cells
<b>SPECIFICITY :</b>	Recognises thrombomodulin (TM), a 105 kDa protein
<b>QUANTITY :</b>	0.1 mg / 1.0 ml
<b>FORMAT :</b>	Phosphate-buffered saline. With 1% BSA and 0,09% Sodium Azide*
<b>STORAGE :</b>	Stable at +2-8°C for 12 months For longer storage, freeze aliquots
<b>APPLICATION :</b>	<b>ELISA Detection Antibody</b>

This Antibody can be used as Detection in a human CD141 sandwich Immunoassay to detect human CD141 in combination with human CD141 Capture Antibody (Cat n° 879.720.001). The suggested coating concentration range below should be optimised by each laboratory for each application.

**ELISA: 0,05-0.5µg/ml**

General ELISA protocol is available at [www.diacclone.com](http://www.diacclone.com)

**Related products : CD141 ELISA Set**

Version 4 - 03.16

For research use only

\*Sodium azide is a poisonous and hazardous substance, handle and dispose of properly

#### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

#### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

#### DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod  
- BP 1985

25020 Besançon cedex - France

Tel: +33 (0)3 81 41 38 38

Fax: +33 (0)3 81 41 36 36

[info@diacclone.com](mailto:info@diacclone.com)



[www.diacclone.com](http://www.diacclone.com)