



Diacclone

DIACLONE Immunology Products

IL-2 ELISA Kit Human IL-2 ELISAKit

CATALOGUE N° :	950.010.048, 950.010.096, 950.010.192
SPECIFICITY :	Recognizes both natural and recombinant human IL-2
RANGE :	31.2 pg / ml - 1000 pg / ml
SENSITIVITY :	< 10 pg / ml
INCUBATION :	1h50
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human cytokines Cross reactivity with simian IL-2
KIT SIZE :	48 (Discovery), 96, 192 tests
REFERENCES :	Altokka-Uzun, G. et al., Cephalalgia, 2015; 33:1024-1027 Azadbakht, L. et al., Diabetes Care, 2007; 30(4):967-73 Blanco B. et al., J. Immunol., 2003 ; 171(2):1070-7 Blanco, B. et al., Blood, 2006; 107(9):3575-83 Buenestado, A. et al., J Parenter Enteral Nutr., 2006; 30(4):286-96 Calabria-Linares, C. et al., J. Cell Sci., 2011; 124(5): 820-830 Casati, C. et al., Cancer Res., 2003; 63(15):4507-15 Doniz-Padilla, L. et al., Eur. J. Endocrinol., 2011 ; 165(1): 129-136. Eilertsen, G. O. et al., Lupus, 2011 ; 20(6): 607-613. Ge, N. L. et al., World J Gastroenterol., 2003; 9(10): 2182-5 Giaccone, G. et al., Clin Cancer Res., 2002; 8 (12): 3702-9 Gritzapis, A. D. et al., Br J Cancer, 2003; 88(8): 1292-300 Guo, Z. et al., Biomed Rep., 2016; 5(5): 567-573 Hadley, E. A. et al., Clin Exp Immunol., 2005; 140(1): 101-8 Hohn, H. et al., Clin Exp Immunol., 2003; 131(1): 102-10 Kapsogeorgou E.K. et al., J. Immunol., 2001; 166(5): 3107 - 3113 Kashyap, B. et al., J Trop Pediatr., 2011; fmr093 Kerr J. et al., J. Gen. Virol., 2001; 82(Pt 12): 3011-3019 Liu, C. P. et al., J Leukoc Biol., 2007; 81: 1276 - 1286 Llorente, L. et al., Arthritis Rheum., 2000; 43(8):1790-800 Martin-Cofreces, N. B. et al., J Immunol., 2006; 176(7):4201-7 Martin-Cofreces, N. B. et al., J. Cell Biol., 2008; 182(5): 951-962. Merlo, A. et al., Infect Immun., 2001; 69(10):6022-9 Montcuquet, N. et al., Immunology, 2008; 125(3): 320-30. Nagao, Y. et al., Malar J., 2008; 7: 113. Ralph, C. et al., Clin. Cancer Res., 2010; 16(5): 1662-1672 Robles-Valero, J. et al., Blood, 2010; 115(21): 4174-4184. Rodriguez-Zapata, M. et al., Infect. Immun., 2010; 78:3272-3279 Rybka, J. et al., Clin Exp Med., 2015; 16(4): 493-502. Saifullah, S. Khan, et al., Chem Cent J., 2013; 7:163 Saverino D. et al., J. immunol., 2002; 168(1): 207 - 215 Shaw, D. M. et al., Br J Cancer, 2007; 96(4): 567-74. Sobjanek, M. et al., Postepy Dermatol Alergol., 2016; 33(4): 263-8. Stebbing, J. et al., Ann Oncol., 2003; 14(11):1660-6 Stebbing, J. et al., Clin Exp Immunol., 2004; 138(2): 312-6. Su, C-C et al., Cancer Res., 2004; 64(6): 2148 - 2152 Sun, Z. et al., Cent Eur J Immunol., 2014; 39(2): 216-22 Tang, Y. et al., Am J Trop Med Hyg., 2008; 79(2): 154-158. Tiemessen M.M. et al., Int. Immunol., 2003; 15(12): 1495 - 1504 Trznadel-Grodzka, E. et al., Postepy Hig Med Dosw (Online), 2012; 66: 843-7 Trznadel-Grodzka, E. et al., Postepy Hig Med Dosw (Online), 2012; 66: 848-54 Urso, K. et al., Blood, 2011; 118: 795 - 803 van der Burg, S. H. et al., PNAS, 2007; 104(29): 12087-12092. Woerly G. et al., J. Exp. Med., 1999; 190(4): 487 - 495
KIT CONTENT :	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Calibrated to WHO 86/504

Version 17 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com

