



Diacclone

DIACLONE Immunology Products

IL-6 HS ELISA Kit

Human IL-6 High Sensitivity ELISA Kit

CATALOGUE N° :	950.035.048, 950.035.096, 950.035.192
SPECIFICITY :	Recognizes both natural and recombinant human IL-6
RANGE :	1.56 pg / ml - 50 pg / ml
SENSITIVITY :	< 0.8 pg / ml
INCUBATION :	3h45
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human cytokines No interference with soluble receptor gp80
KIT SIZE :	48 (discovery) 96, 192 tests
REFERENCES :	Borrione, P. et al., Br. J. Sports Med., 2008; 42(11): 894-900. Botella-Carretero, J. I. et al., Eur J Endocrinol., 2005; 153(2): 223-30 Cassidy, E. M. et al., J Psychopharmacol., 2002; 16(3): 230-4. Giannessi, D. et al., Metabolism, 2011; 60(2): 227-33 Giannessi, D. et al., Cell Stress Chaperones, 2007; 12(3): 265-74 Ikonomidis, I. et al., Circulation, 2008; 117(20): 2662-2669. Klapinska, B. et al., Eur J Appl Physiol., 2013 ; 113(11): 2781-93 Krupski, W. et al., Medicine (Baltimore), 2016; 95(7): e2817. Mukhopadhyay, B. et al., Diagn Res 10(2): BC08-10. Pemberton, P. W. et al., BMC Res Notes, 2009; 2: 83 Pingitore, A. et al., J. Clin. Endocrinol. Metab., 2008; 93:1351-1358 Rahaus M. et al., J. Gen. Virol., 2004; 85(Pt 12): 3529 - 3540 Richards, A. A. et al., Am J Hum Biol., 2013 ; 25(4): 457-64 Sanchez-Lazaro, I. J. et al., Heart Int., 2011; 6(2): e17 Waskiewicz, Z. et al., Eur J Appl Physiol., 2012; 112(5): 1679-88. Wehlin L. et al., Eur. J. cardiothorac. Surg., 2003; 25(1): 35 - 42
KIT CONTENT :	Pre-coated 12 strip plates, Biotinylated Secondary Antibody, Amplifier, Standard, Buffers, Streptavidin-HRP, TMB, Stop Reagent

Calibrated to WHO 89/548

Version 9 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com