



Diacclone

## DIACLONE Immunology Products

### IL-12 p70 ELISA Kit

Human IL-12p70 ELISA Kit

<b>CATALOGUE N° :</b>	950.070.048, 950.070.096, 950.070.192
<b>SPECIFICITY :</b>	Recognizes both natural and recombinant human IL-12 Specific to p70 heterodimer
<b>RANGE :</b>	6.25 pg / ml - 200 pg / ml
<b>SENSITIVITY :</b>	2.2 pg / ml
<b>INCUBATION :</b>	3h40
<b>SAMPLE TYPES :</b>	Serum Plasma Cell culture supernatant
<b>SAMPLE SIZE :</b>	100 µl
<b>CROSS REACTION :</b>	No cross reactivity with other human cytokines
<b>KIT SIZE :</b>	48 (Discovery), 96, 192 tests
<b>REFERENCES :</b>	Belderbos, M E. et al., Clin Immunol.,2009; 133(2): 228-37. Bellanger, A-P. et al., Clin. Vaccine Immunol., 2013 ; 20(8): 1133-1142. Boissel, N. et al., J Leukoc Biol.,2006; 79(4): 747-56. Botella-Carretero J. I. et al., Eur. J. Endocrinol., 2005; 153(2): 223 - 230 Charbonnier, A S. et al., J Leukoc Biol.,2003; 73(1): 91-9. El-Emshtay, H. M. et al., Dis Markers 2015: 707254 Foucher, E. D. et al., PLoS One,2013;8(2):e56045 Gironella J. et al., Gut, 2005; 54(9): 1244 - 1253 Gosset P. et al., J. Immunol., 2003 ; 170(10):4943 - 4952 Guo, Z. et al., Biomed Rep.,2016;5(5): 567-573 Hatzfeld-Charbonnier, A. S. et al., J Leukoc Biol.,2007;81: 1179 - 1187 Hernandez, C. et al., PLoS One,2014 ; 9(6): e98458 Hessle, C. et al., Infect Immun.,2000;68(6): 3581-6. Hessle, C. et al., Clin Exp Immunol.,1999; 116(2): 276-82 Krispin, A. et al., Blood,2006; 108(10): 3580-9. Lakomy, D. et al., J. Immunol., 2011; 187: 2775 - 2782 Lu, C. et al., Lupus, 2015;24(1):18-24 Madej, M. et al., Reumatologia,2015; 53(1): 9-13. Olivar, R. et al., J. Immunol.,2016;196(10): 4274-4290 Op den Brouw, M. L. et al., Immunology,2009 ; 126(2): 280-9 Pichavant, M. et al., J Immunol., 2006;177(9): 5912-9. Ratajczak, C. et al., J Biomed Biotechnol., 2007(1): 71921 Rea D. et al., Blood, 2000; 95(10): 3162 - 3167 Rea D. et al., J. Immunol., 2001; 166(8): 5236 -5244 Rea D. et al., J. Virol., 1999; 73(12): 10245 -10253 Schmitt, C. et al., J Leukoc Biol.,2000; 68(6): 836-44. Sene, D. et al., PLoS Pathog.,2010; 6(11): e1001184 Sotiropoulou, P. A. et al., Stem Cells, 2006; 24(1): 74-85. Stebbing J. et al., Ann. Onc., 2003; 14(11): 1660 - 1666 Stebbing, J. et al., Clin Exp Immunol.,2004; 138(2): 312-6. Surbatovic, M. et al., Sci Rep.,2015;5:11355 Taront, S. et al., Part Fibre Toxicol.,2009;6:9 Verdijk R M. et al., J. Immunol., 1999; 163(1): 57 - 61
<b>KIT CONTENT :</b>	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Calibrated to WHO 95/544

Version 11 - 06.17

For research use only

### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

#### DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod  
- BP 1985  
25020 Besançon cedex - France  
Tel: +33 (0)3 81 41 38 38  
Fax: +33 (0)3 81 41 36 36  
info@diacclone.com



[www.diacclone.com](http://www.diacclone.com)