



Diacclone

## DIACLONE Immunology Products

### IL-12 p40/p70 ELISA Kit

Human IL-12 p40 + p70 ELISA Kit

<b>CATALOGUE N° :</b>	950.075.048, 950.075.096, 950.075.192
<b>SPECIFICITY :</b>	Recognizes the p40 monomer, dimer or heterodimer (p40 + p35) forms of the molecule
<b>RANGE :</b>	62.5 pg / ml - 2000 pg / ml
<b>SENSITIVITY :</b>	< 20 pg / ml
<b>INCUBATION :</b>	1h45
<b>SAMPLE TYPES :</b>	Serum Plasma Cell culture supernatant
<b>SAMPLE SIZE :</b>	100 µl
<b>CROSS REACTION :</b>	No cross reactivity with other human cytokines Cross reacts with feline and canine IL-12 p40
<b>KIT SIZE :</b>	48 (discovery), 96, 192 tests
<b>REFERENCES :</b>	Beyth S. et al., Blood, 2005; 105(5): 2214 - 2219 Bogels, M. et al., Oncoimmunology, 2012; 1(6): 798-809. Braitch, M. et al., Arch Neurol., 2008; 65(5): 633-635. Guerra-Laso, J. M. et al., Immunology, 2015; 144(2): 291-301 Hoonakker, M. E. et al., PLoS One, 2016; 11(8): e0161428. Nagao, Y. et al., Malar J., 2008; 7: 113. Rodriguez-Zapata, M. et al., Infect. Immun., 2010; 78:3272-3279 Sun, Z. et al., J Biomed Res., 2012; 26(6): 456-66 Tenca C. et al., J. Immunol., 2005; 174(11): 6757 - 6763 Vilisaar, J. et al., J Neuroimmune Pharmacol., 2015; 10(3): 457-67
<b>KIT CONTENT :</b>	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Calibrated to WHO 95/544

Version 9 - 06.17

For research use only

### 25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

### Fast Track Your Research...

#### DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAPlex

#### DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod  
- BP 1985  
25020 Besançon cedex - France  
Tel: +33 (0)3 81 41 38 38  
Fax: +33 (0)3 81 41 36 36  
info@diacclone.com



[www.diacclone.com](http://www.diacclone.com)