



Diacclone

DIACLONE Immunology Products

CD25 ELISA Kit

Human CD25 ELISA kit

CATALOGUE N° :	950.500.048, 950.500.096, 950.500.192
SPECIFICITY :	Recognizes both natural and recombinant human soluble IL-2R
RANGE :	68.75 pg / ml - 2220 pg / ml
SENSITIVITY :	< 32.5 pg / ml
INCUBATION :	3h40
SAMPLE TYPES :	Serum Plasma Cell culture supernatant
SAMPLE SIZE :	100 µl
CROSS REACTION :	No cross reactivity with other human soluble receptors
KIT SIZE :	48 (Discovery), 96, 192 test
REFERENCES :	Beijes, M. G. H. et al., Nephrol. Dial. Transplant, 2015; gfv092 Botella-Carretero J. I. et al., Eur. J. Endocrinol., 2005; 153(2): 223 - 230 Christiansson, L. et al., PLoS One, 2013; 8(1): e55818 Dennert, R. et al., Clin. Vaccine Immunol., 2012; 19(8): 1182-1187 Foussat A. et al., Blood, 2001; 98(6): 1678 - 1686 Lindqvist, C. A. et al., Immunology, 2010; 131(3): 371-6 Magnusson, M. K. et al., United European Gastroenterology Journal, 2013: 2050640613502962 Malmstrom, P.-U. et al., Clin. Cancer Res., 2010; 16(12): 3279-3287 Marcellis, R. G. J. et al., Eur. Respir. J., 2011; 38: 628 - 634 Mostard, R. L. et al., BMC Pulm Med., 2012; 12(1): 57 Rasool, R. et al., Int J Organ Transplant Med., 2015; 6(1): 8-13 Sobjanek, M. et al., Postepy Dermatol Alergol., 2016; 33(4): 263-8. Verwoerd, A. et al., Eur. Respir. J., 2015; 09031936.00005315 Wijnen, P. A. et al., Eur. Respir. J., 2014 ; 43(6): 1730-1739 Wilde, B. et al., Arthritis Res Ther., 2016; 18: 147. Zedan, M. et al., ISRN Pediatr., 2013; 2013: 824781 Zekri, A. R. et al., World J Gastroenterol., 2005; 11(42): 6624-30
KIT CONTENT :	Pre-coated 12 strip plates, biotinylated secondary antibody, standards, controls (when available), buffers, Streptavidin-HRP, TMB, Stop Reagent.

Version 15 - 06.17

For research use only

25 YEARS OF EXPERTISE

With over 25 years experience and extensive expertise we are committed to providing excellence in Monoclonal Antibody and Immunoassay development.

The expanding range of Diacclone Immunology products is specifically designed to advance research applications.

Our experience and expertise coupled to the diversity of our product range makes Diacclone a clear choice to

Fast Track Your Research...

DIACLONE PRODUCTS



Monoclonal Antibodies



ELISA



ELISpot



DIAplex

DIACLONE SAS

6 rue Dr Jean-François-Xavier Girod
- BP 1985
25020 Besançon cedex - France
Tel: +33 (0)3 81 41 38 38
Fax: +33 (0)3 81 41 36 36
info@diacclone.com



www.diacclone.com