

Protein Molecular Weight Marker (14 - 120kDa) - MB-201-0200
Code: MB-201-0200

Size: 200 µL


Product Description: Protein Molecular Weight Marker (14 - 120kDa) - MB-201-0200

Concentration: 1X

PhysicalState: Liquid (in 1x Loading Buffer)

Label	Unconjugated
Buffer	See application note.
Stabilizer	None
Preservative	0.01% (w/v) Sodium Azide
Storage Condition	Store vial at -20° C. Expiration date is six (6) months from date of opening container.
Synonyms	Protein ladder, Protein marker, Protein standard, Protein molecular weight standard, Molecular weight marker, Protein molecular weight ladder, Protein marker, Mw marker, Mw standard
Application Note	Protein Molecular Weight Marker is designed for use in accurate estimation of protein molecular weight in SDS-polyacrylamide gel electrophoresis and western blotting. Thaw the cocktail at room temperature. Gently invert vial several times to ensure that the solution is homogeneous and that any precipitated material is re-dissolved. A loading volume of 5-10µl is appropriate for most gel formats.
Background	This protein molecular weight marker is a cocktail containing seven purified proteins that range in size from 14 to 120 kDa and thus span a molecular weight range common to most proteins and protein subunits. The protein mixture has been formulated to yield well-defined bands after SDS-polyacrylamide gel electrophoresis and staining with Coomassie blue or other common staining techniques.
Purity And Specificity	Protein Molecular Weight Marker is supplied in 1X loading buffer and is ready-to-use. It contains: [62.5 mM Tris-HCl, pH 6.8, 2% (w/v) SDS, 250 mM NaCl, 1mM EDTA, 4% (v/v) -Mercaptoethanol, 35% (v/v) Glycerol and 0.05% (w/v) bromophenol blue.] This product are unstained markers optimized for use with Coomassie blue staining, but the loading volume can be adjusted for compatibility with other staining methods (i.e. silver staining). The intensity of the 67 kDa bovine albumin band has been increased to serve as a reference.
Assay Dilutions	User Optimized
Other Assays	User Optimized
Related Products	
	200-301-268 Anti-AKT pS473 (MOUSE) Monoclonal Antibody - 200-301-268
	610-4302 Anti-MOUSE IgG (H&L) (RABBIT) Antibody Peroxidase Conjugated - 610-4302
	611-1302 Anti-RABBIT IgG (H&L) (GOAT) Antibody Peroxidase Conjugated - 611-1302
	B304 NORMAL GOAT SERUM (NGS) - B304


Images



4-20% Tris-Glycine SDS-PAGE
Stained with Coomassie Blue

2

SDS-Page of Protein Molecular Weight Marker. Lane 1: Protein Molecular Weight Marker 12.5 uL. Lane 2: Protein Molecular Weight Marker 10.0 uL. Lane 3: Protein Molecular Weight Marker 7.5 uL. Lane 4: Protein Molecular Weight Marker 5.0 uL. Exact molecular weights are detailed on image 1.



Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.