

## BPP Bioportide™ Bacterial Transformation Kit for *cyanobacteria*

### Product Overview

The BPP Bioportide™ kits provide a ready-to-use system for efficient transformation of various microorganisms. Each vial contains **lyophilized BPP Bioportide™**, supplemented by a separate vial with **sterile 50% glycerol** for storage at -20°C. This innovative reagent facilitates the uptake of nucleic acids without requiring inducing cell stress, i.e. heat shock, or competent cells, i.e. by means of electroporation or chemical treatment of cell wall and cell membrane.

### Components Supplied

1. **BPP Bioportide™**: Contains HEPES and arginine as buffer material. One vial contains **1 µg** of lyophilized BPP Bioportide™ which is equal to **ten reactions** for transformation of *cyanobacteria*. One reaction corresponds to a volume of transformed cells of 1 ml, which is available for plating. Kits are available in sizes 30, 50 and 100 reactions.
2. **BPP Bioportide™ Buffer**: The fluid component is delivered in a separate vial that contains 500 µl of sterile 50% glycerol.

### Storage, Handling and Stability

- Store lyophilized BPP Bioportides™ at **-20°C** or -80°C for long-term storage, respectively.
- The **BPP Bioportide™ buffer** can also be stored at **-20°C** or -80°C.
- Once dissolved, the **BPP Bioportide™ working solution** must be stored at **-20°C**.

### Caution

BPP Bioportides™ should be handled with care. As effects to prolonged or unprotected exposure are not foreseeable, precautions are strictly advised, and users should always follow mandatory lab safety regulations and especially adhere to the following safety standards:

- Avoid skin or eye contact. Wearing gloves and eye protection when handling the product is strongly recommended. In case of contact, rinse skin or eyes thoroughly with water immediately.
- Do not inhale or consume. Avoid the formation of and exposure to aerosols, e.g. by working under a safety cabinet or in well-ventilated workspaces. Avoid ingestion.
- Do not use it while pregnant or breastfeeding.
- Only use within the parameters detailed in this manual and do not leave open on a bench.

## Product Description

**BPP Bioportide™** is a recombinant protein with a molecular weight of roughly **37.5 kDa**. The exact size depends on the kind of BPP Bioportide™. They are distinguished by their cell-penetrating domain, which were designed for specific cell envelope compositions. The protein enables the transport of nucleic acids, ranging from **22 nt to >10 kb**, into cells of a wide range of organisms, i.e. bacteria, archaea, fungi, plants, and animals (depending on the BPP Bioportide™ in use). Both **RNA** and **DNA**, including single- and double-stranded, linear or circular forms, can be transferred reliably. BPP Bioportides™ are designed for robust performance in extreme environments (up to 80°C, pH > 9.0, high salt concentrations).

## Key Features

- No need for competent cells, heat shock, electroporation or any other pretreatment.
- Works with **nanogram** quantities of nucleic acids in **cyanobacteria** like **Synechocystis or Nostoc**.
- Efficient transformation, with up to ~300 CFU for Synechocystis transformed with 10 ng of a 6 kbp linear integrative plasmid (DNA).
- Short and easy protocol. Less time per experiment, less DNA and less cells are needed compared to standard natural transformation.

## General Considerations

- You can adapt our standard protocol as a starting point to the specific needs of your strain. If the BPP Bioportide™ protocol does not offer enough details, proceed according to your standard protocol, as this contains the optimum conditions for your strain. If you are not sure, please do not hesitate to contact us.
- Note, that you do not have to conduct pretreatment of your cells. BPP Bioportide™ transformation works best with early to mid log phase cells. **The optimal OD<sub>750</sub> is 0.4.**
- As the carrier domain has an affinity for polysaccharides, avoid those within the medium that you use for incubation to bring in contact your nucleic acid / BPP Bioportide™ complex with your cells.
- Consider using a known and easy to handle methylation adapted plasmid as control in addition to your main experiment.

## General Transformation Protocol for *cyanobacteria*

### 1. Material preparation

- Grow cells according to your standard protocol until the absorption (750 nm) reached 0.4. Centrifuge 50 ml of cells (3200xg for 10 minutes) and resuspend in 2 ml of appropriate medium (e.g. BG 11).
- At the end, one reaction of BPP Bioportides™ corresponds to a volume of transformed cells of 1 ml, which is available for plating. Please calculate the material for your planned experiments accordingly.
- Prepare the **appropriate medium** (e.g. BG 11) and **agar plates** with and without appropriate antibiotics.

### 2. Prepare BPP Bioportide™ Working Solution

- Add **10 µl of BPP Bioportide™ buffer** to one vial of BPP Bioportide™ (yielding a **100 ng/µl working solution**).
- Incubate at **room temperature** for **20 minutes** without agitating.

### 3. Prepare Nucleic Acid Mixture

- Mix **1 µl** of nucleic acid (10 ng/µl) with **1 µl of BPP Bioportide™ Working solution**.
- Incubate the mixture at **room temperature** for **10 minutes**.

### 4. Transform the Cells

- Use **50 µl of cyanobacterial cell culture** (see step 1).
- Add **50 µl** of cells to the BPP Bioportide™/nucleic acid mixture (**2 µl**; see step 3) and keep at **room temperature** for **10 minutes**.

### 5. Recover and Plate

- Add **950 µl of appropriate medium (e.g. BG 11)** to the cells and keep in continuous shaking and light (~80 µmol) for **60 minutes**. In addition, test the standard incubation time of your strain, e.g. 6 hours.
- Transfer cells to **HATF membrane on BG 11 agar plates without antibiotics** for **24 hours**.
- Transfer membrane to BG 11 agar plates **with appropriate antibiotics**. Colonies appear after approximately **10 days** of incubation.

## Documentation and support

### Customer and technical support

Visit [www.badische-peptide-proteine.de](http://www.badische-peptide-proteine.de) for the latest service and support information.

### Limited product warranty

Badische Peptide und Proteine GmbH warrant their products as set forth in the General Terms and Conditions of Sale at [www.badische-peptide-proteine.de](http://www.badische-peptide-proteine.de). If you have any questions, please contact Badische Peptide und Proteine GmbH at [www.badische-peptide-proteine.de](http://www.badische-peptide-proteine.de).

### Manufacturer

Badische Peptide und Proteine GmbH, Waldshuter Straße 5, 68239 Mannheim, Germany

### Liability

TO THE EXTENT ALLOWED BY LAW, BPP GmbH WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

### Revision

The information in this user guide is subject to change without notice.

Latest version: June 2025

### Licensing and Trademark

IMPORTANT INSTRUCTIONS – READ CAREFULLY: Limited Use Label License (LULL)

This product is covered by U.S., European and other pending patent applications and issued patents. The purchase of this product conveys to the buyer a non-transferable right to use the purchased product for research purposes only. This product may not be resold or incorporated in any manner into another product for resale. Any use for therapeutics or diagnostics is strictly prohibited. Chemical analysis and/or reverse engineering of the product or its components is strictly prohibited. Any use except for research purposes requires a license from Badische Peptide und Proteine GmbH. For further information and license terms and conditions, please contact Badische Peptide und Proteine GmbH (visit [www.badische-peptide-proteine.de](http://www.badische-peptide-proteine.de)).

©2024 Badische Peptide und Proteine GmbH. All rights reserved. All trademarks are the property of Badische Peptide und Proteine GmbH and its subsidiaries unless otherwise specified.